


Assay Class: High Sensitivity DNA Assay
Data Path: C:\...gh Sensitivity DNA Assay_DE24802100_2017-03-13_14-44-45.xad

Created: 3/13/2017 2:44:45 PM
Modified: 4/10/2017 2:15:20 PM

Electrophoresis File Run Summary



Instrument Information:












Instrument Name: DE24802100
Serial#: DE24802100
Firmware: C.01.069
Type: G2938B

Assay Information:

Assay Origin Path: C:\Program Files (x86)\Agilent\2100 bioanalyzer\2100 expert\assays\dsDNA\High Sensitivity DNA.xsy
Assay Class: High Sensitivity DNA Assay
Version: 1.03
Assay Comments: Copyright © 2003-2010 Agilent Technologies

Chip Information:

Chip Lot #:
Reagent Kit Lot #:
Chip Comments:



Assay Class: High Sensitivity DNA Assay
Data Path: C:\...gh Sensitivity DNA Assay_DE24802100_2017-03-13_14-44-45.xad

Created: 3/13/2017 2:44:45 PM
Modified: 4/10/2017 2:15:20 PM

Electrophoresis File Run Summary (Chip Summary)

Sample Name	Sample Comment	Rest. Digest	Status Observation	Result Label	Result Color
1-Gcarol (1:40)	C. Delwiche	<input type="checkbox"/>	✓		
2-Gbelize (1:40)	C. Delwiche	<input type="checkbox"/>	✓		
3-GexcentExp (1:40)	C. Delwiche	<input type="checkbox"/>	✓		
4-GexcentLog (1:40)	C. Delwiche	<input type="checkbox"/>	✓		
sample 5		<input type="checkbox"/>	✓		
sample 6		<input type="checkbox"/>	✓		
1-Gcarol (1:40)		<input type="checkbox"/>	✓		
sample 8		<input type="checkbox"/>	✓		
sample 9		<input type="checkbox"/>	✓		
2-Gbelize (1:40)		<input type="checkbox"/>	✓		
3-GexcentExp (1:40)		<input type="checkbox"/>	✓		
Ladder		<input type="checkbox"/>	✓		

Chip Lot #**Reagent Kit Lot #****Chip Comments :**

Assay Class: High Sensitivity DNA Assay
Data Path: C:\...gh Sensitivity DNA Assay_DE24802100_2017-03-13_14-44-45.xad

Created: 3/13/2017 2:44:45 PM
Modified: 4/10/2017 2:15:20 PM

Electrophoresis Assay Details

General Analysis Settings

Number of Available Sample and Ladder Wells (Max.) : 12
Minimum Visible Range [s] : 32
Maximum Visible Range [s] : 138
Start Analysis Time Range [s] : 33
End Analysis Time Range [s] : 137.5
Ladder Concentration [pg/μl] : 1950
Uses Standard Area for Ladder Fragments
Lower Marker Concentration [pg/μl] : 125
Upper Marker Concentration [pg/μl] : 75
Used Upper Marker for Quantitation
Standard Curve Fit is Point to Point
Show Data Aligned to Lower and Upper Marker

Integrator Settings

Integration Start Time [s] : 33.05
Integration End Time [s] : 137
Slope Threshold : 0.8
Height Threshold [FU] : 5
Area Threshold : 0.1
Width Threshold [s] : 0.6
Baseline Plateau [s] : 0.5

Filter Settings

Filter Width [s] : 0.5
Polynomial Order : 4


Ladder

Ladder Peak	Size	Area
1	35	160
2	50	210
3	100	208
4	150	221
5	200	242
6	300	270
7	400	305
8	500	306
9	600	336
10	700	321
11	1000	366
12	2000	413
13	3000	411
14	7000	400
15	10380	214

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...gh Sensitivity DNA Assay_DE24802100_2017-03-13_14-44-45.xad

Created: 3/13/2017 2:44:45 PM
 Modified: 4/10/2017 2:15:20 PM

Electropherogram Summary



Overall Results for Ladder

Noise: 0.5


Peak table for Ladder

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	50	150.00	4,545.5	Ladder Peak
3	100	150.00	2,272.7	Ladder Peak
4	150	150.00	1,515.2	Ladder Peak
5	200	150.00	1,136.4	Ladder Peak
6	300	150.00	757.6	Ladder Peak
7	400	150.00	568.2	Ladder Peak
8	500	150.00	454.5	Ladder Peak
9	600	150.00	378.8	Ladder Peak
10	700	150.00	324.7	Ladder Peak
11	1,000	150.00	227.3	Ladder Peak
12	2,000	150.00	113.6	Ladder Peak
13	3,000	150.00	75.8	Ladder Peak
14	7,000	150.00	32.5	Ladder Peak
15	10,380	75.00	10.9	Upper Marker

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...gh Sensitivity DNA Assay_DE24802100_2017-03-13_14-44-45.xad

Created: 3/13/2017 2:44:45 PM
 Modified: 4/10/2017 2:15:20 PM

Electropherogram Summary Continued ...



Overall Results for sample 1 : 1-Gcarol (1:40)

Number of peaks found: 17 Corr. Area 1: 736.3
 Noise: 0.5

Peak table for sample 1 : 1-Gcarol (1:40)

Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	242	130.90	820.0	
3	246	38.98	240.0	
4	252	299.24	1,796.2	
5	334	13.74	62.4	
6	350	12.66	54.7	
7	361	24.25	101.9	
8	388	9.81	38.3	
9	429	7.60	26.8	
10	438	6.02	20.8	
11	628	223.01	538.4	
12	694	163.48	356.9	
13	1,128	32.15	43.2	
14	1,454	9.93	10.3	
15	1,663	7.02	6.4	
16	1,774	13.20	11.3	
17	2,072	22.30	16.3	
18	10,380	75.00	10.9	Upper Marker
19	12,904	0.00	0.0	


Region table for sample 1 : 1-Gcarol (1:40)

From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
200	502	736.3	47	311	25.0	617.52	3,235.4	Blue

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...gh Sensitivity DNA Assay_DE24802100_2017-03-13_14-44-45.xad

Created: 3/13/2017 2:44:45 PM
 Modified: 4/10/2017 2:15:20 PM

Electropherogram Summary Continued ...



Overall Results for sample 2 : 2-Gbelize (1:40)

Number of peaks found: 14 Corr. Area 1: 735.3
 Noise: 0.4

Peak table for sample 2 : 2-Gbelize (1:40)


Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	244	448.69	2,785.3	
3	337	23.23	104.5	
4	358	9.19	38.9	
5	378	7.59	30.5	
6	400	6.67	25.3	
7	412	6.21	22.8	
8	442	10.75	36.9	
9	601	170.29	429.2	
10	651	58.37	135.8	
11	680	196.77	438.5	
12	1,601	6.79	6.4	
13	1,789	18.32	15.5	
14	2,411	14.53	9.1	
15	10,380	75.00	10.9	Upper Marker
16	12,964	0.00	0.0	

Region table for sample 2 : 2-Gbelize (1:40)

From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
200	503	735.3	47	312	25.6	581.28	3,049.7	Blue

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...gh Sensitivity DNA Assay_DE24802100_2017-03-13_14-44-45.xad

Created: 3/13/2017 2:44:45 PM
 Modified: 4/10/2017 2:15:20 PM

Electropherogram Summary Continued ...**Overall Results for sample 3 : 3-GexcentExp (1:40)**

Number of peaks found: 10 Corr. Area 1: 715.9
 Noise: 0.3

Peak table for sample 3 : 3-GexcentExp (1:40)


Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	241	164.23	1,031.1	
3	258	217.01	1,275.3	
4	307	67.95	334.8	
5	345	16.60	73.0	
6	406	11.60	43.3	
7	450	6.58	22.1	
8	641	218.61	516.5	
9	653	35.52	82.4	
10	725	170.99	357.2	
11	10,380	75.00	10.9	Upper Marker
12	12,934	0.00	0.0	

Region table for sample 3 : 3-GexcentExp (1:40)

From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
200	498	715.9	49	309	24.4	712.37	3,749.4	Blue

Assay Class: High Sensitivity DNA Assay
 Data Path: C:\...gh Sensitivity DNA Assay_DE24802100_2017-03-13_14-44-45.xad

Created: 3/13/2017 2:44:45 PM
 Modified: 4/10/2017 2:15:20 PM

Electropherogram Summary Continued ...**Overall Results for sample 4 : 4-GexcentLog (1:40)**

Number of peaks found: 11 Corr. Area 1: 759.3
 Noise: 0.3

Peak table for sample 4 : 4-GexcentLog (1:40)


Peak	Size [bp]	Conc. [pg/μl]	Molarity [pmol/l]	Observations
1	35	125.00	5,411.3	Lower Marker
2	50	7.34	222.0	
3	241	148.87	934.3	
4	246	298.48	1,836.4	
5	346	10.25	44.9	
6	654	127.88	296.4	
7	703	32.00	68.9	
8	788	49.02	94.3	
9	1,326	9.41	10.8	
10	1,489	9.09	9.2	
11	1,701	4.78	4.3	
12	10,380	75.00	10.9	Upper Marker
13	12,792	0.00	0.0	

Region table for sample 4 : 4-GexcentLog (1:40)

From [bp]	To [bp]	Corr. Area	% of Total	Average Size [bp]	Size distribution in CV [%]	Conc. [pg/μl]	Molarity [pmol/l]	Color
200	501	759.3	52	308	24.1	741.39	3,902.2	Blue

Assay Class: High Sensitivity DNA Assay
Data Path: C:\...gh Sensitivity DNA Assay_DE24802100_2017-03-13_14-44-45.xad

Created: 3/13/2017 2:44:45 PM
Modified: 4/10/2017 2:15:20 PM

Gel Image

Assay Class: High Sensitivity DNA Assay
Data Path: C:\...gh Sensitivity DNA Assay_DE24802100_2017-03-13_14-44-45.xad

Created: 3/13/2017 2:44:45 PM
Modified: 4/10/2017 2:15:20 PM

Curves

Standard Curve

